



# MicroSquisher

## Micro-Scale Compression System



The MicroSquisher is designed to perform compression testing on specimens between 50 $\mu$ m and 2mm with force resolutions as small as 50nN. Forces are determined from the deflection of a flexible cantilever beam to which one compression plate is attached. Displacement control is achieved by manipulating the base of that beam using a motorized piezo stage.



The specimen can be tested in ambient air or in a temperature-controlled fluid bath. An integrated camera system allows synchronized imaging at up to 5Hz.

## Micro-Scale Material Characterization

The MicroSquisher can be used to determine the compressive stress-strain properties of a variety of materials (hydrogel microspheres, small tissue samples, scaffolds and cell aggregates) with peak forces ranging from 1 $\mu$ N to 1mN.

The device can perform displacement-controlled compression and stress relaxation testing.

*The graphs below demonstrate the range of force measurements that can be achieved using the device. Use different cantilever beam lengths and cross sections to change the force range.*






# MicroSquisher


## Micro-Scale Compression System

### Cell Aggregate Compression

The interface tensions that exist play an important role in the organization of cells within aggregates. These properties can be determined by analyzing the force-time curve and test images from a parallel plate compression test.




The MicroSquisher image analysis module quantifies the aggregate profile, allowing cell-cell and cell-medium interface tensions to be calculated.




### SPECIFICATIONS

Specimen Size ..... 50µm–2mm  
Force Resolution ..... 50nN  
Max. Compression Rate ..... 30µm/s

### System Overview



1. Camera (1280X960 pixel USB)
2. Zoom Lens System (0.3–4.2mm field of view)
3. Vertical Position Stage (manual imaging adjustment)
4. Test Chamber (details below)



1. Temperature controller fluid chamber
2. Vertical stage position feedback sensor (0.1 µm resolution)
3. Axis piezo motor positioning system (0.05µm incremental motion)
4. Cantilever beam load applicator with attached upper compression plate

## World Precision Instruments, Inc.

USA: International Trade Center, 175 Sarasota Center Boulevard, Sarasota FL 34240-9258 USA  
Tel: 941-371-1003 • Fax: 941-377-5428 • E-mail: info@wpiinc.com • Internet: http://www.wpiinc.com

Germany: Zossener Str. 55, 10961 Berlin, Germany • Tel: 030-6188845 • Fax: 030-6188670 • E-mail: wptide@wpi-europe.com

China & Hong Kong: WPI Shanghai Trading Co., Ltd. • Tel: +86 688 85517 • E-mail: ChinaSales@china.wpiinc.com

UK: Astonbury Farm Business Centre • Aston, Stevenage, Hertfordshire SG2 7EG England • Tel: 01438-880025 • E-mail: wpiuk@wpi-europe.com