

Adapter #15978

Warranty

WPI (World Precision Instruments, Inc.) warrants to the original purchaser that this equipment, including its components and parts, shall be free from defects in material and workmanship for a period of one year* from the date of receipt. WPI's obligation under this warranty shall be limited to repair or replacement, at WPI's option, of the equipment or defective components or parts upon receipt thereof f.o.b. WPI, Sarasota, Florida U.S.A. Return of a repaired instrument shall be f.o.b. Sarasota.

The above warranty is contingent upon normal usage and does not cover products which have been modified without WPI's approval or which have been subjected to unusual physical or electrical stress or on which the original identification marks have been removed or altered. The above warranty will not apply if adjustment, repair or parts replacement is required because of accident, neglect, misuse, failure of electric power, air conditioning, humidity control, or causes other than normal and ordinary usage.

To the extent that any of its equipment is furnished by a manufacturer other than WPI, the foregoing warranty shall be applicable only to the extent of the warranty furnished by such other manufacturer. This warranty will not apply to appearance terms, such as knobs, handles, dials or the like.

WPI makes no warranty of any kind, express or implied or statutory, including without limitation any warranties of merchantability and/or fitness for a particular purpose. WPI shall not be liable for any damages, whether direct, indirect, special or consequential arising from a failure of this product to operate in the manner desired by the user. WPI shall not be liable for any damage to data or property that may be caused directly or indirectly by use of this product.

Claims and Returns

Inspect all shipments upon receipt. Missing cartons or obvious damage to cartons should be noted on the delivery receipt before signing. Concealed loss or damage should be reported at once to the carrier and an inspection requested. All claims for shortage or damage must be made within ten (10) days after receipt of shipment. Claims for lost shipments must be made within thirty (30) days of receipt of invoice or other notification of shipment. Please save damaged or pilfered cartons until claim is settled. In some instances, photographic documentation may be required. Some items are time-sensitive; WPI assumes no extended warranty or any liability for use beyond the date specified on the container

Do not return any goods to us without obtaining prior approval and instructions from our Returns Department. Goods returned (unauthorized) by collect freight may be refused. Goods accepted for restocking will be exchanged or credited to your WPI account. Goods returned which were ordered by customers in error are subject to a 25% restocking charge. Equipment which was built as a special order cannot be returned.

Repairs

Contact our Customer Service Department for assistance in the repair of apparatus. Do not return goods until instructions have been received. Returned items must be securely packed to prevent further damage in transit. The Customer is responsible for paying shipping expenses, including adequate insurance on all items returned for repairs. Identification of the item(s) by model number, name, as well as complete description of the difficulties experienced should be written on the repair purchase order and on a tag attached to the item.

* Electrodes, batteries and other consumable parts are warranted for 30 days only from the date on which the customer receives these items.



Adapter #15978

Allows pulling of P-5 and P-7 multi-barrel glass on micropipette puller PUL-100

INSTRUCTION MANUAL

Serial No. _____

070601



World Precision Instruments, Inc.

USA

International Trade Center, 175 Sarasota Center Blvd., Sarasota FL 34240-9258
Tel: 941-371-1003 • Fax: 941-377-5428 • E-mail: sales@wpiinc.com

UK

Astonbury Farm Business Centre • Aston, Stevenage, Hertfordshire SG2 7EG
Tel: 01438-880025 • Fax: 01438-880026 • E-mail: wpiuk@wpi-europe.com

Germany

Liegnitzer Str. 15, D-10999 Berlin
Tel: 030-6188845 • Fax: 030-6188670 • E-mail: wpide@wpi-europe.com

Japan

6-11, Uchikanda 2-Chome, Chiyoda-Ku, Tokyo 101-0047
Tel: 81-3-3258-1641 • Fax: 81-3-3258-1657 • E-mail: info@wpi-j.com

Australia

P.O. Box 1191, Glen Waverley, Victoria 3150
Tel: (03) 9887-6262 • Fax: (03) 9887-9585 • E-mail: wpi.au@ozemail.com.au

Internet

www.wpiinc.com • www.wpi-medical.com
www.nitricoxide.net • www.pipetter.com

World Precision Instruments

FOR PUL-100 MODELS WITH SERIAL
NUMBERS BELOW 42480





Figure 1: PUL-100

Figure 2:
Components of
#15978 Adapter Kit

INSTALLATION

1. Remove both glass clamps ⑥ from the PUL-100.
2. Using the pin press (a modified C-clamp) included in the #15978 Adapter Kit, press pin ⑤ into hole ④ in upper guide block ①.
3. Using the pin press, press pin ⑤ into hole ④ in lower guide block ②.
4. Mount upper block ① onto upper block ① and secure with hold down screws ④.
5. Mount one of the glass clamps ⑥ to upper block.
6. Remove mounting screw from the other glass clamp and mount the spacer ⑥ on it.
7. Mount grab plate ③ to lower block ② and secure to lower block ② with screw and spacer from Step 6.
8. Store all leftover parts for future use in original puller configuration.
9. Mount 3 mm Square Box Filament (WPI #14074) on the puller.



FOR PUL-100 MODELS WITH SERIAL
NUMBERS 42480 AND ABOVE




Figure 1: PUL-100

Figure 2:
Components of
#15978 Adapter Kit

INSTALLATION

1. Remove both glass clamps ⑥ from the PUL-100.
2. Mount upper block ① onto upper block ① and secure with Hold Down Screws ④.
3. Mount one of the glass clamps ⑥ to Upper Block.
4. Remove Mounting Screw from the other glass clamp and mount the spacer ⑥ on it.
5. Mount Grab Plate ③ to Lower Block ② and secure to lower block ② with screw and spacer from Step 4.
6. Bag all leftover parts and place with unit.
7. Mount 3 mm Square Box Filament (WPI #14074) on the puller.

