



WORLD
PRECISION
INSTRUMENTS
Instrumenting scientific ideas

INSTRUCTION MANUAL

MityFlex

Peristaltic Pump

Serial No. _____

www.wpiinc.com

051316

CONTENTS

ABOUT THIS MANUAL 1

INTRODUCTION 1

SPECIFICATIONS..... 2

OPERATING INSTRUCTIONS..... 3

 Installation 3

 Roller/Tubing Selection..... 3

 Variable Speed..... 3

 Tubing Installation 4

 Alternate Tubing Installation Procedure..... 4

 Fuse Installation 5

 Motor Brush Replacement 5

DECLARATION OF CONFORMITY 6

WARRANTY 7

 Claims and Returns 7

 Repairs 7

Copyright © 2016 by World Precision Instruments, Inc. All rights reserved. No part of this publication may be reproduced or translated into any language, in any form, without prior written permission of World Precision Instruments, Inc.

ABOUT THIS MANUAL

The following symbols are used in this guide:



This symbol indicates a CAUTION. Cautions warn against actions that can cause damage to equipment. Please read these carefully.



This symbol indicates a WARNING. Warnings alert you to actions that can cause personal injury or pose a physical threat. Please read these carefully.

NOTES and TIPS contain helpful information.



INTRODUCTION

MITYFLEX™ variable speed drive is powered by DC gearmotors. A specially designed electronic speed control system provides smooth acceleration control over the entire flow range. The three position on/off switch allows the user to maintain speed setting when unit is turned off or reversed. This unit is designed for continuous operation at any speed.

SPECIFICATIONS

ELECTRICAL: 115 VAC, 60 Hz

PRIMING: Unit is self priming and will hold vacuum when off.

OPERATION: Pump can run dry without damaging the motor or drive unit, however the tubing life will be reduced.

CONSTRUCTION: High impact plastic and long lasting alloy metals housed in a steel enclosure.

MAX. HEAD PRESSURE: 40 PSI. Normal pressure 18 PSI varies depending on tubing used.

MAX. SUCTION LIFT: 29 feet (8.8m) with H₂O.

RATING: 100% continuous operation at maximum of 40°C ambient.

FLOW RATE: Variable speed drive system - reversible.

FUSE: Located in the inlet receptacle. Uses 5 X 20mm - 1.0 amp. (See page 4 for replacement procedure.)

DC MOTOR BRUSHES: Replaceable (See page 4 for replacement procedure.)

WEIGHT: 13 lb

Included with MityFlex are four short lengths of Norprene™ tubing to fit into the roller mechanism — $\frac{1}{16}$ " , $\frac{1}{8}$ " , $\frac{3}{16}$ " , $\frac{1}{4}$ " — and black rollers for use with $\frac{3}{16}$ " to $\frac{1}{4}$ " tubing. The optional **MITY-KIT** includes replacement tubing and a set of red rollers for use with $\frac{1}{16}$ " and $\frac{1}{8}$ " tubing.

FLOW RANGE:

Tubing ID	Flow (mL/min)
$\frac{1}{16}$ " (0.062"/1.6mm)	1.7 - 48.0
$\frac{1}{8}$ " (0.125"/3.1mm)	7.0 - 190.0
$\frac{3}{16}$ " (0.187"/4.8mm)	16.0 - 456.0
$\frac{1}{4}$ " (0.250"/6.4mm)	28.0 - 800.0

Flow range based on pumping water at 25°C and 0 PSI using 60 durometer thermoplastic elastomer tubing.

OPERATING INSTRUCTIONS

Installation

- POWER REQUIREMENTS:** Voltage and frequency of power supply must be the same as shown on unit.
- WIRING CONNECTIONS:** All wiring and electrical connections must comply with local and national electrical codes.
- MOUNTING:** Place on any flat surface. Ensure that space is provided at the back of the case for air circulation through ventilation holes. Do not allow tubing to be crimped.
- MOUNTING LOCATION:** Pump should be placed in a dry location with adequate cooling provided. The ambient temperature should not exceed 40°C. ***This unit should not be used outdoors or in hazardous locations.***

Roller/Tubing Selection

Use the following color coded rollers with the appropriate tubing I.D. Your unit is assembled with black rollers. Red rollers are included in the optional **MITY-KIT**. See separate instructions for roller replacement enclosed in your tube kit.

Roller	Tubing ID	Tubing OD	mL/Revolution
Black	$\frac{1}{4}$ " (0.250)	$\frac{7}{16}$ " (0.438"/11.1mm)	3.5
Black	$\frac{3}{16}$ " (0.187)	$\frac{3}{8}$ " (0.375"/9.5mm)	2.1
Red	$\frac{1}{8}$ " (0.125) *	$\frac{3}{8}$ " (0.375"/9.5mm)	0.8
Red	$\frac{1}{16}$ " (0.062)*	$\frac{3}{8}$ " (0.375"/9.5mm)	0.21

- * Red rollers for use with $\frac{1}{8}$ " and $\frac{1}{16}$ " tubing are included in the optional MITY-KIT.
- NOTE:** The $\frac{1}{16}$ " ID tubing is actually a small piece of tubing inside a larger piece.

Variable Speed

- Variable speed control may be set at any level before turning unit on.
- Unit is designed for continuous operation at any speed.








CAUTION: Avoid prolonged operation in the stalled position to ensure longer motor life.

Tubing Installation

Inspect all tubing regularly and replace if any sign of deterioration occurs. Always wear safety glasses and protective clothing when working with chemicals.


1. Disconnect power to the pump.
2. Disconnect the suction and discharge tubing from the pump tubing.
3. Remove four thumb screws and pump cover.
4. Pull out old pump tubing and discard.
5. Clean roller race, removing any particles that could damage tubing.



6. Using extra bracket (Fig. 1), rotate roller bracket assembly until it is in position as in Fig. 2. (Be sure power is disconnected.)
7. Push new tubing into inlet port, anchoring the tubing in grippers (see Fig. 3). If using 9" piece of tube, $1\frac{1}{2}$ " should be extended outside of inlet port.
8. Continue to rotate the roller bracket assembly using the extra bracket, while pushing the tubing into the roller race (Fig. 4).
9. Insert the tubing into the outer port (Fig. 5) and replace the cover and screws.

Alternate Tubing Installation Procedure

Without removing the pump cover, run the unit at high RPM. Pull the old tubing out from outlet port of pump. Install new tube by inserting into inlet port; continue feeding tube $\frac{1}{4}$ " in at a time allowing the roller assembly to push tubing through. When tubing is half way around the race, put your finger over the inlet hole of the tube and continue to feed tubing $\frac{1}{4}$ " at a time; this should make installation easier. Continue until desired tubing extension is obtained. Run unit one minute.



Fuse Installation

1. Disconnect power to pump.
2. Find the power socket on the back of the case.
3. Pull the power cord out.
4. Slide the drawer below the receptacle out toward you.
5. Remove old fuse and replace with new fuse.
6. Push drawer in until it locks in place.

Motor Brush Replacement

1. Disconnect power to pump.
2. Remove case cover by removing 8 screws.
3. With a flat head screw driver, remove the brush caps on both sides of the DC motor.
4. Remove old brush and replace with new brush. (Slide in horizontally)
5. Push metal plate on the end of the brush inside. Hold it while replacing brush cap.
6. Replace case cover.

DECLARATION OF CONFORMITY



WORLD PRECISION INSTRUMENTS, INC.

175 Sarasota Center Boulevard
Sarasota, FL 34240-9258 USA
Telephone: (941) 371-1003 Fax: (941) 377-5428
e-mail wpi@wpiinc.com

DECLARATION OF CONFORMITY

We: World Precision Instruments, Inc.
175 Sarasota Center Boulevard
Sarasota FL 34240-9258
USA

as the distributor of the apparatus listed, declare that the product:

Title: 913 Series MITYFLEX Variable Speed Peristaltic Pump

to which this declaration relates is in conformity with the following standards or other normative documents:


EN61326:2006
EN55011:1998+A2:2002
EN61000-3-2:2001
EN61000-3-3:2001
EN61000-6-2:2001

and therefore conform(s) with the protection requirements of Council Directive 2004/108/EC relating to electromagnetic compatibility.

Issued on: November 6, 2008



Mr. Cliff Bredenberg
General Manager
World Precision Instruments, Inc.
175 Sarasota Center Boulevard
Sarasota, FL 34243-9258 USA



Mr. Glen Carlquist
Vice President of Manufacturing
World Precision Instruments, Inc.
175 Sarasota Center Boulevard
Sarasota, FL 34243-9258 USA

WARRANTY

WPI (World Precision Instruments, Inc.) warrants to the original purchaser that this equipment, including its components and parts, shall be free from defects in material and workmanship for a period of one year* from the date of receipt. WPI's obligation under this warranty shall be limited to repair or replacement, at WPI's option, of the equipment or defective components or parts upon receipt thereof f.o.b. WPI, Sarasota, Florida U.S.A. Return of a repaired instrument shall be f.o.b. Sarasota.

The above warranty is contingent upon normal usage and does not cover products which have been modified without WPI's approval or which have been subjected to unusual physical or electrical stress or on which the original identification marks have been removed or altered. The above warranty will not apply if adjustment, repair or parts replacement is required because of accident, neglect, misuse, failure of electric power, air conditioning, humidity control, or causes other than normal and ordinary usage.

To the extent that any of its equipment is furnished by a manufacturer other than WPI, the foregoing warranty shall be applicable only to the extent of the warranty furnished by such other manufacturer. This warranty will not apply to appearance terms, such as knobs, handles, dials or the like.

WPI makes no warranty of any kind, express or implied or statutory, including without limitation any warranties of merchantability and/or fitness for a particular purpose. WPI shall not be liable for any damages, whether direct, indirect, special or consequential arising from a failure of this product to operate in the manner desired by the user. WPI shall not be liable for any damage to data or property that may be caused directly or indirectly by use of this product.

Claims and Returns

Inspect all shipments upon receipt. Missing cartons or obvious damage to cartons should be noted on the delivery receipt before signing. Concealed loss or damage should be reported at once to the carrier and an inspection requested. All claims for shortage or damage must be made within ten (10) days after receipt of shipment. Claims for lost shipments must be made within thirty (30) days of receipt of invoice or other notification of shipment. Please save damaged or pilfered cartons until claim is settled. In some instances, photographic documentation may be required. Some items are time-sensitive; WPI assumes no extended warranty or any liability for use beyond the date specified on the container

Do not return any goods to us without obtaining prior approval and instructions from our Returns Department. Goods returned (unauthorized) by collect freight may be refused. Goods accepted for restocking will be exchanged or credited to your WPI account. Goods returned which were ordered by customers in error are subject to a 25% restocking charge. Equipment which was built as a special order cannot be returned.

Repairs

Contact our Customer Service Department for assistance in the repair of apparatus. Do not return goods until instructions have been received. Returned items must be securely packed to prevent further damage in transit. The Customer is responsible for paying shipping expenses, including adequate insurance on all items returned for repairs. Identification of the item(s) by model number, name, as well as complete description of the difficulties experienced should be written on the repair purchase order and on a tag attached to the item.

** Electrodes, batteries and other consumable parts are warranted for 30 days only from the date on which the customer receives these items.*

USA

International Trade Center, 175 Sarasota Center Blvd., Sarasota FL 34240-9258
Tel: 941-371-1003 • Fax: 941-377-5428 • E-mail: sales@wpiinc.com

UK

1 Hunting Gate, Hitchin, Hertfordshire SG4 0TJ
Tel: 44 (0)1462 424700 • Fax: 44 (0)1462 424701 • E-mail: wpiuk@wpi-europe.com

Germany

Zossener Str. 55, 10961 Berlin
Tel: 030-6188845 • Fax: 030-6188670 • E-mail: wpide@wpi-europe.com

China & Hong Kong

WPI Shanghai Trading Co., Ltd.
Rm 29a, No8 Dongfang Rd., Pudong District, Shanghai, 200120 PR China
Tel: +86 21 6888 5517 • E-mail: chinasales@china.wpiinc.com

Brazil

Av. Conselheiro Nébias, 756 sala 2611, Santos-CEP: 11045-002, São Paulo Brazil •
Tel: (013) 406-29703 • E-mail: info@brazil.wpiinc.com

Internet

www.wpiinc.com • www.wpi-europe.com • www.wpiinc.cn