



PolyFil

*For coupling multi-barrel microinjection
pipettes to the injection system*

Claims and Returns

• Inspect all shipments upon receipt. Missing cartons or obvious damage to cartons should be noted on the delivery receipt before signing. Concealed loss or damage should be reported at once to the carrier and an inspection requested. All claims for shortage or damage must be made within 10 days after receipt of shipment. Claims for lost shipments must be made within 30 days of invoice or other notification of shipment.

Please save damaged or pilfered cartons until claim settles. In some instances, photographic documentation may be required. Some items are time sensitive; WPI assumes no extended warranty or any liability for use beyond the date specified on the container.

• WPI cannot be held responsible for items damaged in shipment en route to us. Please enclose merchandise in its original shipping container to avoid damage from handling. We recommend that you insure merchandise when shipping. The customer is responsible for paying shipping expenses including adequate insurance on all items returned.

• Do not return any goods to WPI without obtaining prior approval and instructions (RMA#) from our returns department. Goods returned unauthorized or by collect freight may be refused. The RMA# must be clearly displayed on the outside of the box, or the package will not be accepted. Please contact the RMA department for a request form.

• Goods returned for repair must be reasonably clean and free of hazardous materials.

• A handling fee is charged for goods returned for exchange or credit. This fee may add up to 25% of the sale price depending on the condition of the item. Goods ordered in error are also subject to the handling fee.

• Equipment which was built as a special order cannot be returned.

• Always refer to the RMA# when contacting WPI to obtain a status of your returned item.

• For any other issues regarding a claim or return, please contact the RMA department

Warning: This equipment is not designed or intended for use on humans.



World Precision Instruments, Inc.

USA

International Trade Center, 175 Sarasota Center Blvd., Sarasota FL 34240-9258
Tel: 941-371-1003 • Fax: 941-377-5428 • E-mail: sales@wpiinc.com

UK


Astonbury Farm Business Centre • Aston, Stevenage, Hertfordshire SG2 7EG
Tel: 01438-880025 • Fax: 01438-880026 • E-mail: wpiuk@wpi-europe.com

Germany

Liegnitzer Str. 15, D-10999 Berlin
Tel: 030-6188845 • Fax: 030-6188670 • E-mail: wpide@wpi-europe.com


Internet

www.wpiinc.com • www.wpi-medical.com
www.nitricoxide.net • www.pipetter.com



World Precision Instruments

The PolyFil kit includes a plastic multi-barrel micropipette holder/handle, seven **MF28G MicroFil** needles, seven sections of tubing with male luer lock fittings, one flow-through manifold with five luer lock ports, one three-way stopcock, one hot melt glue gun and three glue sticks.



Procedure for making connections:

1. First pull the multi-barrel pipette with a micropipette puller. Mark each barrel with a different color so that they can be identified later.

2. Load each barrel with the solution to be injected. To prevent the liquid from overheating when bonding

the multi-barrel with the hot melt adhesive and prevent liquid crossover between each barrel, make sure that at least 15-20 mm of the top section of multi-barrel pipette is left unfilled.

3. Insert a MicroFil needle into each barrel about 7-10 mm deep, as shown in Fig. 1.

4. Install a stick of hot melt adhesive into the hot melt glue gun. Plug the glue gun into an electric outlet. Let it warm up for five minutes or until the glue can be dispensed. Then bond the MicroFil to the multi-barrel micropipette as shown in Fig. 1.

5. Mount the micropipette holder to a micromanipulator and mount the multi-barrel to the micropipette holder

6. If a seven-barrel micropipette is used, attach the three-way stopcock to one end of the five-port manifold so that the seven barrels can be controlled independently.

7. Connect the tubing as show in Fig. 2.


8. Turn on only the channels of the manifold from which the liquid will be injected. Make sure all the channels that do not need injection are in the off position.

The system is ready to run.

After use, the MicroFil needles can be recovered by immersing the bonded joint in turpentine for ten minutes. Turpentine may be purchased at any paint or hardware store.

The MicroFil syringe needles can be purchased separately from WPI under the part number **MF28G-30**.

Fig. 2



Warranty

WPI (World Precision Instruments, Inc.) warrants to the original purchaser that this equipment, including its components and parts, shall be free from defects in material and workmanship for a period of one year* from the date of receipt. WPI's obligation under this warranty shall be limited to repair or replacement, at WPI's option, of the equipment or defective components or parts upon receipt thereof f.o.b. WPI, Sarasota, Florida U.S.A. Return of a repaired instrument shall be f.o.b. Sarasota.

The above warranty is contingent upon normal usage and does not cover products which have been modified without WPI's approval or which have been subjected to unusual physical or electrical stress or on which the original identification marks have been removed or altered. The above warranty will not apply if adjustment, repair or parts replacement is required because of accident, neglect, misuse, failure of electric power, air conditioning, humidity control, or causes other than normal and ordinary usage.

To the extent that any of its equipment is furnished by a manufacturer other than WPI, the foregoing warranty shall be applicable only to the extent of the warranty furnished by such other manufacturer. This warranty will not apply to appearance terms, such as knobs, handles, dials or the like.

WPI makes no warranty of any kind, express or implied or statutory, including without limitation any warranties of merchantability and/or fitness for a particular purpose. WPI shall not be liable for any damages, whether direct, indirect, special or consequential arising from a failure of this product to operate in the manner desired by the user. WPI shall not be liable for any damage to data or property that may be caused directly or indirectly by use of this product.

* Electrodes, batteries and other consumable parts are warranted for 30 days only from the date on which the customer receives these items.