



MPM-7

Multi-Parameter Monitoring Device



- Monitor seven core parameters, including ECG, heart rate, NIBP, pulse oximetry, respiration rate, apnea trends and temperature
- High resolution, colorful TFT display

MPM-7 Multi-Parameter Monitor offers portable, cost-effective, versatile and reliable multi-parameter monitoring in the animal market. Seven core parameters and options allow the monitor to be used standalone in small animal clinics and large veterinary facilities by connecting, wire or wireless, with a central monitoring system.

The high-resolution color screen clearly displays time, alarms, ECG, heart rate, NIBP, pulse oximetry, respiration rate, apnea trends, temperatures and trace freeze simultaneously. The operators can ceaselessly browse the data in

the previous 240 hours and store 72 hour of data in memory. With optional internal thermal printer, 3 waveforms can be recorded.

To cope with modern clinics equipped with many high-powered electronic instruments, the monitor is designed for protection from defibrillator and high-frequency electrosurgical generator when used on the same animal at the same time. In the event of an unexpected power shortage, a built-in rechargeable battery keeps the monitor continuously working for 2 hours and can be automatically and rapidly recharged.

World Precision Instruments, Inc.

USA: International Trade Center, 175 Sarasota Center Boulevard, Sarasota FL 34240-9258 USA
Tel: 941-371-1003 • Fax: 941-377-5428 • E-mail: info@wpiinc.com • Internet: <http://www.wpiinc.com>

Germany: Zossener Str. 55, 10961 Berlin, Germany • Tel: 030-6188845 • Fax: 030-6188670 • E-mail: wptide@wpi-europe.com

China & Hong Kong: WPI Shanghai Trading Co., Ltd. • Tel: +86 688 85517 • E-mail: ChinaSales@china.wpiinc.com

UK: Astonbury Farm Business Centre • Aston, Stevenage, Hertfordshire SG2 7EG England • Tel: 01438-880025 • E-mail: wpiuk@wpi-europe.com