



WORLD
PRECISION
INSTRUMENTS
Instrumenting scientific ideas

Premixslip® Lubricant & Rust Inhibitor

A ten minute soak keeps surgical instruments looking new



Features

- High-grade mineral oil lubricant and rust inhibitor
- Steam penetrable, steam sterilizable and ETO sterilizable
- Non-oily, non-sticky, non-toxic and silicone free.
- 16 oz. Spray bottle

Benefits

- Prevents rust, spots, stains and corrosion, prolonging the life of the instrument
- Prevents damage due to misuse of stiff or grinding instruments
- Routine use enhances the life and efficiency of instruments
- Greatly reduces repair and replacement cost

Applications

- Lubricate and inhibit rust on all your surgical instruments

Premixslip is a non-sticky, non-toxic, silicone free and water soluble mineral oil lubricant that can withstand the high temperatures of autoclaving without becoming gummy. Forming a protective barrier on surgical instruments, Premixslip prevents rusting, staining and spotting while restoring articulation to sticky box-locks and joints.

WPI Part Number	Description
WP34552	Premixslip Instrument Lubricant and Rust Inhibitor, 16 oz. spray bottle

WORLD PRECISION INSTRUMENTS

USA: International Trade Center, 175 Sarasota Center Boulevard, Sarasota FL 34240-9258 USA
Tel: (941) 371-1003 • Fax: (941) 377-5428 • E-mail: wpi@wpiinc.com • Internet: www.wpiinc.com

UK: 1 Hunting Gate, Hitchin, Hertfordshire SG4 0TJ England • Tel: 44 (0)1462 424700 • E-mail: wpiuk@wpi-europe.com

Germany: Saarstraße 23, D-61169 Friedberg (Hesson) • Tel: +49 (0)6031 1602171 • Fax: +49 (0)6031 1602180 • E-mail: wpide@wpi-europe.com

China & Hong Kong: Rm 25e, No8 Donfang Rd., Pudong District, Shanghai 200120 PRC • Tel: +86 688 85517 • E-mail: ChinaSales@china.wpiinc.com

Brazil: Conselheiro Nabias, 756 sala2611, Santos-São Paulo 11045-002 Brazil • E-mail: info@brazil.wpiinc.com