


EVOM™ MANUAL ELECTRODES

WHICH ELECTRODE IS RIGHT FOR MY APPLICATION?




STX4

GENERAL PURPOSE
TEER ELECTRODE
SKU: EVM-EL-03-03-01



STX HTS

HIGH THROUGHPUT
SCREENING PROBES
SKU: EVM-EL-03-02-XX



ENDOHM

CELL CULTURE
CUP CHAMBERS
SKU: EVM-EL-03-01-XX

IDEAL FOR HIGH
THROUGHPUT SCREENING
(HTS) PLATES



MEASURE INSIDE WELL PLATE



REMOVABLE BLADES



MOST ACCURATE IN
CAPTURING SMALLER CHANGES



HANDS-FREE MEASUREMENT



COMPATIBLE WITH OLDER
EVOMS (REQUIRES ADAPTER)



96-WELL PLATE



24-WELL PLATE



6 & 12-WELL PLATE



**Accurate &
Cost-effective**

**High Throughput
Screening**

**Most Accurate
1% Tolerance**