


INSPECTION CHECKLIST



Individual Parts Labeled



SI-MB4 Parts List		
Base Unit	(1) Base plate	A
	(4) Lower support rods	B
	(4) Upper support rods	C
	(1) Upper horizontal stabilizer bar	D
	(1) Lower horizontal stabilizer bar	E
Reservoirs	(4) Buffer reservoir assemblies include:	F
	Buffer reservoir	
	Reservoir clamp	
	Tubing	
	(Units are labeled 1-4. Unit 1 has aeration manifold built in the clamp and a water flow meter attached. Tubing length varies on each unit depending on its position in they system.)	
Tissue Chambers	(4) Tissue chamber assemblies include:	G
	Tissue chamber	
	Tissue chamber clamp	
	Tissue holder with oxygenation frit	
	Tubing	
Startup Kit	(Units are labeled 1-4. The length of tubing attached to each assembly varies depending on the position of the assembly. The clamp on chamber 4 has two aeration outlets on the right side.)	
	(4) Transducer positioners	H
	(1) SI-OVO2F Vascular tissue holder kit (pkg of 10)	I
	(4) M6 screws	J
	(4) M8 screws with washers	J
	(1) Tool kit which includes:	K
	18mm wrench	
	10mm wrench	
	3 hex wrenches	
	(1) Silicone overflow drain system	L
	(4) Stimulation cables	M
	(4) SI-BAM21-LC transducer amplifiers	N
	(4) SI-KG-20 Force transducers	O
	(1) Instruction Manual	

