



# Hotplates, Magnetic Stirrer



- Ceramic work surface, 7.5 x 7.5 inches
- Stir speed from 60–1500RPM
- Temperature control up to 380°C
- Two models: heat-stir, heat only

This hotplate and hotplate with a magnetic stirrer feature a durable, chemical resistant, white, ceramic work surface. These units easily fit inside a bio-hood or on a crowded counter top. Precise speed and temperature controls offer quick adjustments. 8 x 9" footprint



The magnetic stirrer (top) and the hotplate stirrer (bottom) both feature a ceramic work surface.

## MAGNETIC STIRRER SPECIFICATIONS

Speed Range .....	60–1500RPM
Platform Size .....	7.5 x 7.5" (19 x 19cm)
Control .....	Quick Adjustment Knob
Dimensions .....	8 x 9 x 4.5" (20 x 23 x 11.5cm) WxDxH
Electrical .....	115V or 230V, 50–60 Hz, 300W
Warranty .....	2 years

## HOTPLATE/MAGNETIC STIRRER SPECIFICATIONS

Speed Range .....	60–1500RPM
Temperature Range .....	50–380°C
Platform Size .....	7.5 x 7.5" (19 x 19cm)
Control .....	Quick Adjustment Knob
Dimensions .....	8 x 9 x 4.5" (20 x 23 x 11.5cm) WxDxH
Electrical .....	115V or 230V, 50–60 Hz, 500W
Warranty .....	2 years



## World Precision Instruments, Inc.

USA: International Trade Center, 175 Sarasota Center Boulevard, Sarasota FL 34240-9258 USA

Tel: 941-371-1003 • Fax: 941-377-5428 • E-mail: info@wpiinc.com • Internet: http://www.wpiinc.com

Germany: Zossener Str. 55, 10961 Berlin, Germany • Tel: 030-6188845 • Fax: 030-6188670 • E-mail: wptide@wpi-europe.com

China & Hong Kong: WPI Shanghai Trading Co., Ltd. • Tel: +86 688 85517 • E-mail: ChinaSales@china.wpiinc.com

UK: 1 Hunting Gate, Hitchin, Hertfordshire SG4 0TJ England • Tel: 44 (0)1462 424700 • E-mail: wpiuk@wpi-europe.com



# Hotplates, Magnetic Stirrer



- Ceramic work surface, 7.5 x 7.5 inches
- Stir speed from 60–1500RPM
- Temperature control up to 380°C
- Two models: heat-stir, heat only

This hotplate and hotplate with a magnetic stirrer feature a durable, chemical resistant, white, ceramic work surface. These units easily fit inside a bio-hood or on a crowded counter top. Precise speed and temperature controls offer quick adjustments. 8 x 9" footprint

The magnetic stirrer (top) and the hotplate stirrer (bottom) both feature a ceramic work surface.

## MAGNETIC STIRRER SPECIFICATIONS

Speed Range .....	60–1500RPM
Platform Size.....	7.5 x 7.5" (19 x 19cm)
Control.....	Quick Adjustment Knob
Dimensions .....	8 x 9 x 4.5" (20 x 23 x 11.5cm) WxDxH
Electrical.....	115V or 230V, 50–60 Hz, 300W
Warranty .....	2 years

## HOTPLATE/MAGNETIC STIRRER SPECIFICATIONS

Speed Range.....	60–1500RPM
Temperature Range.....	50–380°C
Platform Size.....	7.5 x 7.5" (19 x 19cm)
Control.....	Quick Adjustment Knob
Dimensions .....	8 x 9 x 4.5" (20 x 23 x 11.5cm) WxDxH
Electrical.....	115V or 230V, 50–60 Hz, 500W
Warranty .....	2 years



## World Precision Instruments, Inc.

**USA:** International Trade Center, 175 Sarasota Center Boulevard, Sarasota FL 34240-9258 USA  
**Tel:** 941-371-1003 • **Fax:** 941-377-5428 • **E-mail:** info@wpiinc.com • **Internet:** http://www.wpiinc.com  
**Germany:** Zossener Str. 55, 10961 Berlin, Germany • Tel: 030-6188845 • Fax: 030-6188670 • E-mail: wpide@wpi-europe.com  
**China & Hong Kong:** WPI Shanghai Trading Co., Ltd. • Tel: +86 688 85517 • E-mail: ChinaSales@china.wpiinc.com  
**UK:** 1 Hunting Gate, Hitchin, Hertfordshire SG4 0TJ England • Tel: 44 (0)1462 424700 • E-mail: wpiuk@wpi-europe.com