

HEMOSTATIC FORCEPS

UNDERSTAND THEIR USES




Bleeding Control

Use to clamp blood vessels and control bleeding during surgery.



Precise Tissue Manipulation

Manipulate delicate tissues, grasp and hold onto structures, and provide stability during procedures.



Atraumatic Clamping

Delicate serrations help prevent tissue damage when clamping blood vessels or other structures.



Hemostasis Assessment

Assess the effectiveness of hemostasis after procedures.



Dissection

Use for dissection to separate tissues. Helpful during microdissection or microsurgery.



Emergency Hemostasis

Valuable in emergency situations for control of bleeding in trauma cases.