



WORLD
PRECISION
INSTRUMENTS

UltraDose® Enzyme Plus Ultrasonic Cleaning Solution

Clean and disinfect in one easy step with dual action UltraDose® Hospital Grade Ultrasonic Cleaner Concentrate

Welcome to the power of three with UltraDose® Enzyme Plus Ultrasonic Cleaning Solution. Designed to work with an Ultrasonic Cleaning System, you'll get broad spectrum cleaning due to our blend of three unique enzymes and surfactants that targets and breaks down proteins, carbohydrates, and fats. Enzyme Plus removes organic debris, residues, biofilm, inorganic soils, tissue, mucous, and blood from instruments, and general medical equipment. Plus, it's safe on stainless steel, metals, glass, plastics, rubber, and aluminum. Enzyme Plus delivers excellent cleaning results for every cleaning cycle.

***Unleash the
cleaning power
of enzymes***

Product Features

- Key part of infection prevention and instrument reprocessing
- Super concentrated, 1 oz. makes on gallon of RTI solution
- Measured, easy-to-use pour dispensing system
- pH balanced and free rinsing formula
- A blend of three unique enzymes that delivers protein-targeted cleaning
- Can also be used as a pre-soak

THE POWER OF 3




32 oz. Bottle
Product Code: WPI038 (U.S. Only)



UltraDose®

ENZYME PLUS

Ultrasonic Cleaning Solution



It's not truly clean until it's UltraDose® clean. UltraDose Enzyme Plus attacks and breaks down protein chains and suspends them for fast and easy ultrasonic cleaning. The result is a thorough cleaning you can count on every time.

Don't trust your instruments with just any cleaner. UltraDose Enzyme Plus attacks protein chains on a molecular level. It's the Enzyme Plus difference. Target problems and get a thorough cleaning—We provide the 360° cleaning solution that professionals rely on when cleaning matters most.

PROTEIN CHAINS


Protein chains including organic debris, blood, and mucous are formed on instruments.

BREAK DOWN

Enzymatic action attacks and breaks down tough protein chains.

SUSPEND AND REMOVE

Proteins are suspended in ultrasonic cleaning solution for easy removal.



Our Ultrasonic Cleaners



Scooba-S Ultrasonic Cleaner Without Heat 4.8l

Order code **WP1-1-ART-A1-00-S**



Digital Ultrasonic Cleaner

Order code **504495**



Economy Ultrasonic Cleaner, 110 V

Order code **UBATH-Y**



Quantrex Ultrasonic Cleaning System

Order code **504216**

WORLD PRECISION INSTRUMENTS