



WORLD  
PRECISION  
INSTRUMENTS

## UltraDose® Germicidal Hospital Grade Ultrasonic Cleaner

Clean and disinfect in one easy step with dual action UltraDose® Hospital Grade Ultrasonic Cleaner Concentrate.

*Don't Take infection  
control for granted!*

Available in 16 oz. containers.

Each container makes  
16 gallons of solution.



The unique blend of an ultrasonic cleaning solution and an E.P.A. registered hospital grade disinfectant come together in this effective product. UltraDose® Germicidal eliminates dangerous microorganisms while simultaneously ultrasonically cleaning and disinfecting instruments, apparatus and equipment. Help maintain the safety of your office and staff and keep your instruments clean and disinfected with UltraDose® Germicidal.

### COVID-19 Effective\*


\*Emerging Viral Pathogen Claim

L&R UltraDose® Germicidal Hospital Grade Ultrasonic Cleaner Concentrate has demonstrated effectiveness against viruses similar to COVID-19 Coronavirus on hard, non-porous surfaces. Therefore, UltraDose® Germicidal Hospital Grade Ultrasonic Cleaner Concentrate can be used against COVID-19 Coronavirus when used in accordance with the directions for use against Feline Calicivirus (surrogate for Norovirus) on hard, non-porous surfaces.

### Product Features

- E.P.A. registered Product
- For use in an ultrasonic cleaner
- Disinfects instruments and the ultrasonic bath
- Cleans non-porous instrument and appliance surfaces, apparatus, and equipment
- Performs as an ultrasonic cleaning solution, disinfectant, deodorant, germicidal, fungicidal, virucidal and tuberculoidal all-in-one
- Concentrated (one ounce makes one gallon.)
- Non-irritating and fragrance-free

## UltraDose® GERMICIDAL HOSPITAL GRADE ULTRASONIC CLEANER



- ✓ Fragrance-Free
- ✓ Super concentrate
- ✓ Dual purpose proprietary formulation
- ✓ Promotes office staff safety

UltraDose® Germicidal is proven effective against viruses much tougher to inactivate than Novel Coronavirus, such as Feline Calicivirus (surrogate for Norovirus). This means UltraDose® Germicidal can be used against COVID-19 Coronavirus on hard, non-porous surfaces.

**Active Ingredient:** Phenylphenol

**Container Size:** 16 oz.

**Dilution:** 1oz./gallon

**Gallons/container:** Makes 16 gallons of ready to use solution

**E.P.A. Registration #:** 34810-31-72869

**Product Code:** WPI036 (U.S. Only)



- COVID-19
- HIV-1
- Pseudomonas Acrusinosa
- Trichophyton
- Salmonella
- Staphylococcus
- Seratia Marcescens
- Herpes Simplex 2
- Mycobacterium Bovis (TB)

### Our Ultrasonic Cleaners



Scooba-S Ultrasonic  
Cleaner Without  
Heat 4.8l

Order code WP1-1-  
ART-A1-00-S



Digital Ultrasonic  
Cleaner

Order code 504495



Economy Ultrasonic  
Cleaner, 110 V

Order code UBATH-Y



Quantrex Ultrasonic  
Cleaning System

Order code 504216

#### WORLD PRECISION INSTRUMENTS

United States: (866) 606-1974 · wpi@wpiinc.com · [www.wpiinc.com](http://www.wpiinc.com)  
UK: +44 (0)1462 424700 · wpiuk@wpi-europe.com · [www.wpi-europe.com](http://www.wpi-europe.com)

Germany: +49 (0)30 61 88845 · wpide@wpi-europe.com · [www.wpi-europe.com](http://www.wpi-europe.com)  
China: +86 21 6888 5517 · chinasales@china.wpiinc.com · [www.wpiinc.net](http://www.wpiinc.net)