



# SURGICAL SCISSORS

## 8 TIPS FOR KEEPING THEM SHARP



### Store Scissors Appropriately

Store your scissors in a clean and dry area. Avoid contact with other instruments. Consider using instrument trays with silicone separators.




### Handle Scissors Mindfully

Avoid dropping them or subjecting them to excessive force.



### Clean Surgical Scissors Promptly

Clean scissors immediately after each use. Wipe the blades with a clean, dry cloth after each use. Use a cloth dampened with a cleaning solution to remove debris.



### Use Scissors Properly

Avoid cutting materials that can dull the blades quickly. Never cut paper, clips, wires, or hard plastics. Avoid cutting materials that are too hard.



### Have Scissors Sharpened Professionally

Have scissors professionally sharpened to ensure clean cuts and quick recoveries.




### Inspect Surgical Scissors Routinely

Inspect scissors for any signs of damage, corrosion, or wear. If you notice any issues, replace them.

### Lubricate Scissors Regularly

Use a surgical instrument lubricant. Work the lubricant in the pivot point of the scissors. Wipe away any excess.



### Sterilize Surgical Scissors Correctly

Sterilize in an open position. Use sterilization baskets or wrap the instruments properly. Follow recommended standards and guidelines.