



# MOUSE ANESTHESIA MASKS

## STEREOTAXIC FRAME COMPATIBILITY CHART


### Mouse Anesthesia Mask

Order Code	Type	Stereotaxic Compatibility	Features	Frame
 <p>505225</p>	Mouse and Rat	505200 Stereotaxic, Warm Base	W mR	ST
		505201 Digital Stereotaxic, Warm Base	W DI mR	
		505206 Dual Stereotaxic, Warm Base	W 2 mR	
		505207 Digital Dual Stereotaxic, Warm Base	W DI 2 mR	
		505205 Mouse & Rat Stereotaxic, Warm Base Only	W mR	
		505203 Ultra Precise Stereotaxic, Warm Base	W mR	
		505204 Ultra Precise Digital Stereotaxic	DI mR	
		505208 Ultra Precise Dual Stereotaxic, Warm Base	W 2 mR	
		505209 Ultra Precise Digital Dual Stereotaxic	DI 2 mR	
		505202 Motorized Stereotaxic, Warm Base	W mR	

### Anesthesia Mask for Mouse Stereotaxic

Order Code	Type	Stereotaxic Compatibility	Features	Frame	
 <p>505226</p>	Mouse and Rat, Neonatal Rat	505281 Mouse & Rat Stereotaxic with Adaptors	mR	ST	
		505283 Mouse & Rat Dual Stereotaxic with Adaptors	2 mR		
		505373 Mouse & Rat Digital Dual Stereotaxic with Adaptors	DI 2 mR		
		505282 Mouse & Rat Ultra Precise Stereotaxic	mR		
		505284 Mouse & Rat Ultra Precise Dual Stereotaxic	mR		
		505233 Mouse & Neonatal Rat Stereotaxic with Adaptors	mn		O
		505381 Mouse & Neonatal Rat Motorized	mn		
		505243 Mouse & Neonatal Rat Adaptor	mn		
	Mouse	505318 Mouse Stereotaxic, Base only	m	ST	
		505313 Mouse Stereotaxic	m		
		505314 Mouse Digital Stereotaxic	DI m		
		505319 Mouse Dual Stereotaxic	2 m		
		505320 Mouse Digital Dual Stereotaxic	DI 2 m		
		505316 Mouse Ultra Precise Stereotaxic	m		UP
		505317 Mouse Ultra Precise Digital Stereotaxic	DI m		
		505321 Mouse Ultra Precise Dual Stereotaxic	2 m		
		505322 Mouse Dual Digital Ultra Precise Stereotaxic	DI 2 m		
		505315 Mouse Motorized Stereotaxic	m		M

### Mouse Stereotaxic Anesthesia Mask with Bite Bar

Order Code	Type	Stereotaxic Compatibility	Features	Frame
 <p>505385</p>	Mouse	505388 Mouse Stereotaxic, Warm Base	W m	ST
		505389 Mouse & Neonatal Rat Digital Stereotaxic with Adaptors	DI mn	
	Mouse and Rat	505200 Stereotaxic, Warm Base	W mR	ST
		505201 Digital Stereotaxic, Warm Base	W DI mR	
		505206 Dual Stereotaxic, Warm Base	W 2 mR	
		505207 Digital Dual Stereotaxic, Warm Base	W DI 2 mR	
		505205 Mouse & Rat Stereotaxic, Warm Base Only	W mR	
		505203 Ultra Precise Stereotaxic, Warm Base	W mR	
		505204 Ultra Precise Digital Stereotaxic	DI mR	
		505208 Ultra Precise Dual Stereotaxic, Warm Base	W 2 mR	
		505209 Ultra Precise Digital Dual Stereotaxic	DI 2 mR	
	505202 Motorized Stereotaxic, Warm Base	W mR	M	
	Rat *with 505242 Adaptor	505211 Standard Stereotaxic Instrument	R	ST
		505216 Dual Stereotaxic	2 R	
		505212 Ultra Precise Stereotaxic	R	
		505217 Ultra Precise Dual Stereotaxic	2 R	
		505293 Motorized Stereotaxic, 18°	R	

ST Standard  
 UP Ultra Precise  
 M Motorized  
 O Other  
 W Warming  
 DI Digital  
 2 Dual Arm  
 R Rat  
 mR Mouse & Rat  
 m Mouse  
 mn Mouse & Neonatal