

# TISSUE VS DRESSING FORCEPS

## UNDERSTANDING THEIR DIFFERENCES

### TYPE



TISSUE FORCEPS



DRESSING FORCEPS

### APPLICATIONS



Use to hold and manipulate delicate tissues during surgical procedures, laboratory work, and dental procedures.

Use to hold and manipulate dressings and other materials during wound care or other medical procedures such as suturing.



Microsurgery



Wound Dressing



Ophthalmic Surgery



Suturing



Dental Surgery



Wound Debridement

### DESIGN



- Fine, serrated jaws for a precise grip on tissues without causing damage.

- Long, thin handles and a wider, flatter tip to prevent damage to dressings and wound sites

- Has the option of 1x2, 2x3, 3x4, or 4x5 teeth.

- May be smooth tipped, without any teeth or serrations, or they may be serrated.



1x2



2x3



3x4



4x5



Sharp



Blunt



Serrated